

READ BEFORE USE

Warning

- Only use parts and accessories provided. Parts and accessories from unauthorised sellers may cause damage and will be not be covered under warranty.
- THIS INCLUDES the charging cable (output 5.0V, 1.0A).
- Only abide by the troubleshooting guide. If the device needs servicing beyond the troubleshooting guide, stop using and contact Bird Healthcare immediately.
- DO NOT attempt to service yourself unless specified.
- DO NOT fill the medication cup beyond 10ml – this can cause electric shock and/or damage to the device.
- DO NOT use without liquid.
- DO NOT use with thick/sticky (viscous) liquid as this can block the mesh. ONLY use with Salbutamol and saline solution.
- Avoid getting the device wet – DO NOT expose the device to running water or submersion.
- DO NOT use the device in hot environments (>40°C).
- DO NOT use the device near flammable gas or near oxygen and anaesthetic mixtures.
- DO NOT expose the device to corrosive liquids and gases.
- DO NOT wrap the power cord around the device.
- DO NOT subject the device and accessories to strong vibrations and impacts.
- DO NOT touch the centre of the spray mesh by hand or with any other objects.

Not recommended for people with:

- Pulmonary Embolism
- Pulmonary Infarction or Odema

If using for **acute asthma** and you feel you are not improving, seek urgent medical advice.

- For the type, dose, and regime of your medication, be sure to follow the instructions from your doctor.
- For children and persons with special needs, use with supervision. Use the provided mask if required.
- This device is to be used for the specified purpose of medication nebulisation. DO NOT use the device for any other purpose.
- Clean and disinfect the medication cup and accessories before using for the first time.
- Keep the medication cup empty and dry when not in use.
- This device is intended for single patient use to avoid cross infection.

Safety notes

DO NOT use without liquid
ONLY use with Salbutamol and/or Saline solution
Replace Medication Cup every 3-6 months

Warranty

e-chamber Portable

Mesh Nebuliser Pro

Product Code: BHC-NEBPRO-1



BIRD HEALTHCARE

For all warranty claims, contact Bird Healthcare directly at sales@birdhealthcare.com or call us at TOLL FREE 1300 365 561 with the following information and a copy of your receipt/invoice

Date of purchase* _____

Customer name* _____

Contact number* _____

Where Purchased _____

Presentation of this card will be required when a warranty service is requested. A new card will not be issued if this card is lost

The warranty period of the body of the device is 12 months.

The warranty period of the medication cup is 3 months.

The masks and mouthpiece are not covered under warranty.

Warranty claims will be reviewed on a case by case basis and we will replace or repair parts if deemed necessary.

Issues not covered by warranty include:

- Damage by human behaviours: touching the mesh & soaking/wetting the device
- Damage caused by not correctly operating the product according to the Quick Start Guide
- Damage caused by accident
- User disassembly without authorisation

Please keep the packaging of this product so that it can be used to return the device when maintenance is required.

Submit to sales@birdhealthcare.com or call us at TOLL FREE 1300 365 561

Bird Healthcare Pty Ltd
 18 Corporate Blvd
 Bayswater VIC Australia 3153



Quick Start Guide



Portable Nebuliser
 with travel case,
 masks & mouthpiece




BIRD HEALTHCARE
 www.birdhealthcare.com

e-chamber Portable Nebuliser Pro

Product Code: BHC-NEBPRO-1

Contents

- a. Lid x 1
- b. Device x 1
- c. Medication cup x 1
- d. Mouthpiece x 1
- e. Mask x 2
- f. Travel case x 1
- g. USB charging cable x 1



Use is for liquid medication nebulisation only

First use

1. Clean, disinfect and dry all of the parts (See *Cleaning*).
2. Charge the device (See *Charging*).

Charging

- Charge the device using the provided USB charging cable – if using a standard adapter, make sure it complies with IEC 60601-1, 5.0V.
- The device is usable whilst charging.
- The device takes approximately 1 hour to fully charge.
- When the device is on low battery it will flash slowly.
- Please charge the device at least once per month, even when not in use.

Cleaning

1. Remove the medication cup, lid, mouthpiece or mask from the device.
2. Submerge and hand wash in luke warm water (no more than 50°C) for about 5 minutes.
3. Set aside to dry in a clean place and if required then disinfect (See *Disinfect*).

Disinfect (if required)

Place cleaned medication cup, lid and mouthpiece in disinfectant solution (such as a Milton per instructions). Or steam sterilise in a standard baby steam device as per device instructions.

Air dry completely after disinfecting

Storage: Store the device in the provided carry case ensuring the following:

- The temperature is between -2 0°C and 55°C.
- The relative humidity is less than 80%.
- Atmospheric pressure is between 70 and 106 kPa.
- The device is not in storage for 5+ years unused.


Instructions For Use

1



Remove the lid on the medication cup.

2




Fill to below the max 10ml capacity line.

3



Close the lid.

4



Attach the mouth piece or mask if desired.

5



Click the blue ON/OFF switch to start (pull the slider down if concealing the switch).

6



Hold the nebuliser vertical to keep the liquid in contact with mesh.

- 7 Take slow deep breathes to inhale the medication.
- 8 After 10 minutes, the device will automatically shut off. If longer use is desired, simply switch the device on again ensuring there is enough liquid medication.
- 9 **NB:** When finished, nebulise with 3ml of warm tap water. Once a week, nebuliser with 3ml of warm tap water and 1ml of vinegar. Press ON/OFF switch to turn off.
- 10 Detach the medication cup by pushing the 'press' button behind the device and pull the medication cup upwards. Dispose of any medication.
- 11 Clean any accessories that were used during nebulisation (See *Cleaning*).

Troubleshooting

Trouble	Cause	Solution
Doesn't turn on and light isn't on	Device isn't on	Try turning on the device and ensure the light comes on
	Button may be jammed	Try pressing ON/OFF switch to dislodge
	Device isn't charging	See <i>Charging</i>
	USB Charging Cable may be broken	See <i>Warranty</i> or use a cable with the same voltage and amp readings (this may void your warranty)
Atomisation is too small	There isn't enough medication	Refill no further than the 10ml line
	The wrong medication is being used	See <i>Warnings</i>
	The spray mesh is blocked	Place 3 drops of white vinegar and 3-6ml of water inside the medication cup, then fully atomise (DO NOT INHALE), then clean the medication cup
There is medication residual	This is normal	
There is a strange noise	There is a fault in the product	See <i>Warranty</i>